

Product Data Sheet

Recombinant Human ATAT1 Protein

Catalog #	Source	Reactivity	Applications
CST1426	E. coli	H	
Description	Recombinant Human ATAT1 Protein is produced by E. coli expression system and the target gene encoding Asp311-Arg410 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ATAT1		
Alternative Names	C6orf134; MEC17; Alpha-tubulin N-acetyltransferase 1; Alpha-TAT; Alpha-TAT1; TAT; Acetyltransferase mec-17 homolog		
Entrez Gene	79969		
SwissProt	Q5SQI0		
Purity			
Chemical Structure	DPGGSPAQRR RTRGTPPGLV AQSCCYSRHG GVNSSSPNTG NQDSKQGEQE TKNRSASEEQ ALSQDGSGEK PMHTAPPQAP APPAQSWTVG GDILNARFIR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com