

Product Data Sheet

Recombinant Human SEC23IP Protein

Catalog #	Source	Reactivity	Applications
CST1428	E. coli	H	
Description	Recombinant Human SEC23IP Protein is produced by E. coli expression system and the target gene encoding Asp91-Met440 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SEC23IP		
Alternative Names	SEC23-interacting protein; p125		
Entrez Gene	11196		
SwissProt	Q9Y6Y8		
Purity			
Chemical Structure	DPFGNIGQSP LTTAATSVGQ SGFPKPLTAL PFTTGSQDVS NAFSPSISKA QPGAPPSSLM GINSYLPSQP SSLPPSYFGN QPQGIPQPGY NPYRHTPGSS RANPYIAPPQ LQQCQTPGPP AHPPPSGPPV QMYQMPPGSL PPVPSSVQSP AQQQVPARPG APSVQVPSPF LLQNQYEPVQ PHWFYCKEVE YKQLWMPFSV FDSLNLLEIY NSVQPDPEV SVLGTDGGRYD VY		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com