

Product Data Sheet

Recombinant Human AGBL2 Protein

Catalog #	Source	Reactivity	Applications
CST1434	E. coli	H	
Description	Recombinant Human AGBL2 Protein is produced by E. coli expression system and the target gene encoding Met671-Gly870 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	AGBL2		
Alternative Names	CCP2; Cytosolic carboxypeptidase 2; ATP/GTP-binding protein-like 2		
Entrez Gene	79841		
SwissProt	Q5U5Z8		
Purity			
Chemical Structure	MKFTQCLAEL KELLRQEIHK KFHELQQDVD LEGSWSDISL SDIESSTSGS DSSLSDGLPV HLANIADELT QKKKMFKKKK KKSQTRKQR NEQYQKKNLN QKLKLTEDTS EKAGFASTLQ KQPTFFKNSE NSSFLPMKNE NPRLNETNLN RRDKDTPLDP SMATLILPKN KGRMQNKKPG FTVSCSPKRT INSSQEPAPG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com