

Recombinant Human HLTf Protein

Catalog #	Source	Reactivity	Applications
CST1436	E. coli	H	
Description	Recombinant Human HLTf Protein is produced by E. coli expression system and the target gene encoding Lys211-Ala350 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HLTf		
Alternative Names	HIP116A; RNF80; SMARCA3; SNF2L3; ZBU1; Helicase-like transcription factor; DNA-binding protein/plasminogen activator inhibitor 1 regulator; HIP116; RING finger protein 80; SWI/SNF-related matrix-associated actin-dependent regulator of chromatin subfamily		
Entrez Gene	6596		
SwissProt	Q14527		
Purity			
Chemical Structure	KTEFDKLFED LKEDDKTHEM EPAEAIETPL LPHQKQALAW MVSRENSKEL PPFWEQRNDL YYNTITNFSE KDRPENHVG ILADDMGLGK TLTAIAVILT NFHDGRPLPI ERVKKNNLLKK EYNVNDDSMK LGGNNTSEKA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com