

Recombinant Human SCRIB Protein

Catalog #	Source	Reactivity	Applications
CST1437	E. coli	H	
Description	Recombinant Human SCRIB Protein is produced by E. coli expression system and the target gene encoding Glu661-Asp1190 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SCRIB		
Alternative Names	CRIB1; KIAA0147; LAP4; SCRB1; VARTUL; Protein scribble homolog; Scribble; hScrib; Protein LAP4		
Entrez Gene			
SwissProt	Q14160		
Purity			
Chemical Structure	EEEEEEGSPQ EEEEEEEEN RAEEEASTE EEDKEGAVVS APSVKGVSFD QANNLLIEPA RIEEEELTLT ILRQTGGLGI SIAGGKGSTP YKGDDEGIFI SRVSEEGPAA RAGVRVGDKL LEVNGVALQG AEHHEAVEAL RGAGTAVQMR VWRERMVEPE NAVTITPLRP EDDYSPRERR GGGRLRPLLP PESPGPLRQR HVACLARSER GLGFSIAGGK GSTPYRAGDA GI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com