

Product Data Sheet

Recombinant Human MEIS2 Protein

Catalog #	Source	Reactivity	Applications
CST1442	E. coli	H	
Description	Recombinant Human MEIS2 Protein is produced by E. coli expression system and the target gene encoding Met61-Lys310 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MEIS2		
Alternative Names	MRG1; Homeobox protein Meis2; Meis1-related protein 1		
Entrez Gene	4212		
SwissProt	O14770		
Purity			
Chemical Structure	MPASMGSAVN DALKRDKDAI YGHPLFPLLA LVFEKCELAT CTPREPGVAG GDVCSSDSFN EDIAVFAKQV RAEKPLFSSN PELDNLMIQA IQVLRFHLL LEKVHELCDN FCHRYISCLK GKMPIDLVID ERDGSSKSDH EELSGSSTNL ADHNPSSWRD HDDATSTHSA GTPGPSSGGH ASQSGDNSSE QGDGLDNSVA SPGTGDDDDP DKDKKRQKKR GIFPKVATNI MR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com