

Product Data Sheet

Recombinant Human SNRPF Protein

Catalog #	Source	Reactivity	Applications
CST1448	E. coli	H	
Description	Recombinant Human SNRPF Protein is produced by E. coli expression system and the target gene encoding Leu11-Glu80 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SNRPF		
Alternative Names	PBSCF; Small nuclear ribonucleoprotein F; snRNP-F; Sm protein F; Sm-F; SmF		
Entrez Gene	6636		
SwissProt	P62306		
Purity			
Chemical Structure	LNGLTGKPVM VKLKWGM EYK GYLVSVDGYM NMQLANTEEY IDGALSGHLG EVLIRCNNVL YIRGV EEEEE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com