

Product Data Sheet

Recombinant Human SARNP Protein

Catalog #	Source	Reactivity	Applications
CST1450	E. coli	H	
Description	Recombinant Human SARNP Protein is produced by E. coli expression system and the target gene encoding Leu11-Lys200 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SARNP		
Alternative Names	HCC1; SAP domain-containing ribonucleoprotein; Cytokine-induced protein of 29 kDa; Nuclear protein Hcc-1; Proliferation-associated cytokine-inducible protein CIP29		
Entrez Gene	84324		
SwissProt	P82979		
Purity			
Chemical Structure	LKLAELKQEC LARGLETGKI KQDLIHLRLQA YLEEHAEAAA NEEDVLGDET EEEETKPIEL PVKEEEPPEK TVDVAAEKKV VKITSEIPQT ERMQKRAERF NVPVSLESKK AARAARFGIS SVPTKGLSSD NKPMVNLDKL KERAQRFGLN VSSISRKSED DEKLKKRKER FGIVTSSAGT GTTEDETEAKK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com