

Product Data Sheet

Recombinant Human SEC24C Protein

Catalog #	Source	Reactivity	Applications
CST1453	E. coli	H	
Description	Recombinant Human SEC24C Protein is produced by E. coli expression system and the target gene encoding Pro501-Arg740 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SEC24C		
Alternative Names	KIAA0079; Protein transport protein Sec24C; SEC24-related protein C		
Entrez Gene	9632		
SwissProt	P53992		
Purity			
Chemical Structure	PPAFIFMIDV SYNAIRTGLV RLLCEELKSL LDFLPREGGA EESAIRVGfV TYNKVLHFYN VKSSLAQPQM MVVSDVADMF VPLLDGFLVN VNESRAVITS LLDQIPEMFA DTRETETVFV PVIQAGMEAL KAAECAGKLF LFHTSLPIAE APGKLKNRDD RKLINTDKEK TLFQPQTGAY QTLAKECVAQ GCCVDLFLFP NQYVDVATLS VVPQLTGGSV YKYASFQVEN DQ		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com