

Product Data Sheet

Recombinant Human ECPAS Protein

Catalog #	Source	Reactivity	Applications
CST1455	E. coli	H	
Description	Recombinant Human ECPAS Protein is produced by E. coli expression system and the target gene encoding Gly1671-Ser1770 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ECM29		
Alternative Names	KIAA0368; Proteasome-associated protein ECM29 homolog; Ecm29		
Entrez Gene	23392		
SwissProt	Q5VYK3		
Purity			
Chemical Structure	GVRTTKNEEE NEKEKELQLE YLLGAFESLG KAWPRNAETQ RCYRQELCKL MCERLKLSTW KVQLGVLQSM NAFFQGLMLL EEEHADPEAL AEILLETCKS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com