

Product Data Sheet

Recombinant Human KIF1C Protein

Catalog #	Source	Reactivity	Applications
CST1464	E. coli	H	
Description	Recombinant Human KIF1C Protein is produced by E. coli expression system and the target gene encoding Thr431-Ser680 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIF1C		
Alternative Names	KIAA0706; Kinesin-like protein KIF1C		
Entrez Gene	10749		
SwissProt	O43896		
Purity			
Chemical Structure	TESQIGPEEA MERLQETEKI IAEINETWEE KLRKTEALRM EREALLAEMG VAVREDGGTV GVFSPKKTPH LVNLNEDPLM SECLLYHIKD GVTRVGQVDM DIKLTGQFIR EQHCLFRSIP QPDGEVVVTL EPCEGAETYV NGKLVTEPLV LKSGNRIVMG KNHVFRFNHP EQARLERERG VPPPPGPPSE PVDWNFAQKE LLEQQGIDIK LEMEKRLQDL ENQYRKEKEE AD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com