

Product Data Sheet

Recombinant Human GNG11 Protein

Catalog #	Source	Reactivity	Applications
CST1465	HEK293	H	
Description	Recombinant Human GNG11 Protein is produced by HEK293 expression system and the target gene encoding Pro2-Ser73 is expressed with N-terminus GST & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GNG11		
Alternative Names	GNGT11; Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-11		
Entrez Gene	2791		
SwissProt	P61952		
Purity			
Chemical Structure	PALHIEDLPE KEKLKMEVEQ LRKEVKLQRQ QVSKCSEEEK NYIEERSGED PLVKGIPEDK NPFKEKGSCV IS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com