

Product Data Sheet

Recombinant Human RPN1 Protein

Catalog #	Source	Reactivity	Applications
CST1471	E. coli	H	
Description	Recombinant Human RPN1 Protein is produced by E. coli expression system and the target gene encoding Asn501-Ser590 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	RPN1		
Alternative Names	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit 1; Dolichyl-diphosphooligosaccharide--protein glycosyltransferase 67 kDa subunit; Ribophorin I; RPN-I; Ribophorin-1		
Entrez Gene	6184		
SwissProt	P04843		
Purity			
Chemical Structure	NRYKQSRDIS TLNSGKKSLE TEHKALTSEI ALLQSRLKTE GSDLCDRVSE MQKLDAQVKE LVLKSAVEAE RLVAGKLKKD TYIENEKLIS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com