

Recombinant Human CYTIP Protein

Catalog #	Source	Reactivity	Applications
CST1472	E. coli	H	
Description	Recombinant Human CYTIP Protein is produced by E. coli expression system and the target gene encoding Glu81-Gly290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CYTIP		
Alternative Names	PSCDBP; Cytohesin-interacting protein; Cytohesin binder and regulator; CYBR; Cytohesin-associated scaffolding protein; CASP; Cytohesin-binding protein HE; Cbp HE; Pleckstrin homology Sec7 and coiled-coil domains-binding protein		
Entrez Gene	9595		
SwissProt	O60759		
Purity			
Chemical Structure	EKQDNETFGF EIQSYRPQNQ NACSSEMFTL ICKIQEDSPA HCAGLQAGDV LANINGVSTE GFTYKQVVDL IRSSGNLLTI ETLNGTMILK RTELEAKLQV LKQTLKQKWV EYRSLQLQEH RLLHGDAANC PSLENMDLDE LSLFGPLPGP GPALVDRNRL SSESSCKSWL SSMTMDSEDG YQTCVSEDSS RGAFSRQTST DDECIFIPKEG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com