

Product Data Sheet

Recombinant Human CSDE1 Protein

Catalog #	Source	Reactivity	Applications
CST1474	E. coli	H	
Description	Recombinant Human CSDE1 Protein is produced by E. coli expression system and the target gene encoding Asn11-Asp290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CSDE1		
Alternative Names	D1S155E; KIAA0885; NRU; UNR; Cold shock domain-containing protein E1; N-ras upstream gene protein; Protein UNR		
Entrez Gene	7812		
SwissProt	O75534		
Purity			
Chemical Structure	NGHNGYPNGT SAALRETGVI EKLTSYGFI QCSERQARLF FHCSQYNGNL QDLKVGDDVE FEVSSDRRTG KPIAVKLVKI KQEILPEERM NGQVVCAVPH NLEKSPAAP GQSPTGSVCY ERNGEVFYLT YTPEDVEGNV QLETGDKINF VIDNNKHTGA VSARNIMLLK KKQARCQGVV CAMKEAFGFI ERGDVVKEIF FHYSEFKGDL ETLQPGDDVE FTIKDRNGKE VA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com