

Product Data Sheet

Recombinant Human CPSF6 Protein

Catalog #	Source	Reactivity	Applications
CST1476	E. coli	H	
Description	Recombinant Human CPSF6 Protein is produced by E. coli expression system and the target gene encoding Glu21-Ala210 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CPSF6		
Alternative Names	CFIM68; Cleavage and polyadenylation specificity factor subunit 6; Cleavage and polyadenylation specificity factor 68 kDa subunit; CFIm68; CPSF 68 kDa subunit; Pre-mRNA cleavage factor Im 68 kDa subunit; Protein HPBR11-4/7		
Entrez Gene	11052		
SwissProt	Q16630		
Purity			
Chemical Structure	EA EYGGHDQI DLYDDVISPS ANNGDAPEDR DYMDTLPTTV GDDVGKGAAP NVVYTYTGKR IALYIGNLTW WTTDEDL TEA VHS LGVNDIL EIKFFENRAN GQSKGFALVG VGSEASSKKL MDLLPKRELH GQNPVVTPCN KQFLSQFEMQ SRKTTQSGQM SGEGKAGPPG GSSRAAFPQG GRGRGRFPGA		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com