

Product Data Sheet

Recombinant Human IKZF3 Protein

Catalog #	Source	Reactivity	Applications
CST1492	E. coli	H	
Description	Recombinant Human IKZF3 Protein is produced by E. coli expression system and the target gene encoding Lys11-Met270 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	IKZF3		
Alternative Names	ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3		
Entrez Gene	22806		
SwissProt	Q9UKT9		
Purity			
Chemical Structure	KSTQEQSVPA ESAAVLNDYS LTKSHEMENV DSGEGPANED EDIGDDSMKV KDEYSERDEN VLKSEPMGNA EEPEIPYSYS REYNEYENIK LERHVVVSFDS SRPTSGKMNC DVCGLSCISF NVLMVHKRSH TGERPFQCNQ CGASFTQKGN LLRHIKLHTG EKPFKCHLCN YACQRRDALY GHLRTHSVEK PYKCEFCGRS YKQRSSLEEY KERCRTFLQS TDPGDTASAE AR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com