

Product Data Sheet

Recombinant Human CHORDC1 Protein

Catalog #	Source	Reactivity	Applications
CST1493	E. coli	H	
Description	Recombinant Human CHORDC1 Protein is produced by E. coli expression system and the target gene encoding Gly11-Glu150 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CHORDC1		
Alternative Names	CHP1; Cysteine and histidine-rich domain-containing protein 1; CHORD domain-containing protein 1; CHORD-containing protein 1; CHP-1; Protein morgana		
Entrez Gene	26973		
SwissProt	Q9UHD1		
Purity			
Chemical Structure	GQRFDPETNS DDACTYHPGV PVFHDALKGW SCCKRRTTDF SDFLSIVGCT KGRHNSEKPP EPVKPEVKTT EKKELCCLKP KFQEHIIQAP KPVEAIKRPS PDEPMTNLEL KISASLKQAL DKLKLSSGNE ENKKEEDNDE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com