

Product Data Sheet

Recombinant Human TCF20 Protein

Catalog #	Source	Reactivity	Applications
CST1495	E. coli	H	
Description	Recombinant Human TCF20 Protein is produced by E. coli expression system and the target gene encoding Ile1261-Leu1780 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TCF20		
Alternative Names	KIAA0292; SPBP; Transcription factor 20; TCF-20; Nuclear factor SPBP; Protein AR1; Stromelysin-1 PDGF-responsive element-binding protein; SPRE-binding protein		
Entrez Gene	6942		
SwissProt	Q9UGU0		
Purity			
Chemical Structure	ISPIPSKRQS QDVKNSSTED KGRLLHSSKE GADKAFNSYA HLSHSQDIKS IPKRDSSKDL PSPDSRNCPA VTLTSPAKTK ILPPRKGRGL KLEAIVQKIT SPNIRRSASS NSAEAGGDTV TLDDILSLKS GPPEGGSVAV QDADIEKRKG EVASDLVSPA NQELHVEKPL PRSSEEWGRS VDDKVKTETH AETVTAGKEP PGAMTSTTSQ KPGSNQGRPD GSLGGTAPLI FP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com