

Product Data Sheet

Recombinant Human EPS15L1 Protein

Catalog #	Source	Reactivity	Applications
CST1496	E. coli	H	
Description	Recombinant Human EPS15L1 Protein is produced by E. coli expression system and the target gene encoding Tyr401-Asp820 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EPS15L1		
Alternative Names	EPS15R; Epidermal growth factor receptor substrate 15-like 1; Eps15-related protein; Eps15R		
Entrez Gene	58513		
SwissProt	Q9UBC2		
Purity			
Chemical Structure	YSLEQDIREK EEAIRQKTSE VQELQNLDLR ETSSLQELEA QKQDAQDRLD EMDQQKAKLR DMLSDVRQKC QDETQMISSL KTQIQSQESD LKSQEDDLNR AKSELNRLQQ EETQLEQSIQ AGRVQLETII KSLKSTQDEI NQARSKLSQL HESRQEAHRS LEQYDQVLDG AHGASLTDLA NLSEGVSLAE RGSFGAMDDP FKNKALLFSN NTQELHPDPF QTEDPFKSDP FK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com