

Product Data Sheet

Recombinant Human UVRAG Protein

Catalog #	Source	Reactivity	Applications
CST1497	E. coli	H	
Description	Recombinant Human UVRAG Protein is produced by E. coli expression system and the target gene encoding Ser391-Phe690 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UVRAG		
Alternative Names	UV radiation resistance-associated gene protein; p63		
Entrez Gene	7405		
SwissProt	Q9P2Y5		
Purity			
Chemical Structure	SRSTIKDNIN DKLTEKEREFLYPLKGGGEKL QFDYGVYLLN KNIAQLRYQH GLGTPDLRQT LPNLKNFMEH GLMVRCDRHH TSSAIPVPKR QSSIFGGADV GFSGGIPSPD KGHRKRASSE NERLQYKTPP PSYNSALAQP VTTVPSMGET ERKITSLSSS LDTSLDFSKE NKKKGEDLVG SLNGGHANVH PSQEQGEALS GHRATVNGTL LPSEQAGSAS VQLPGEFHPV SE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com