

Recombinant Human SLAIN2 Protein

Catalog #	Source	Reactivity	Applications
CST1498	E. coli	H	
Description	Recombinant Human SLAIN2 Protein is produced by E. coli expression system and the target gene encoding Asp241-Arg560 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SLAIN2		
Alternative Names	KIAA1458; SLAIN motif-containing protein 2		
Entrez Gene	57606		
SwissProt	Q9P270		
Purity			
Chemical Structure	DRNPPLSPQS SIDSELSASE LDEDSIGSNY KLNDVTDVQI LARMQEESLR QEYAATTSSRR SSGSSCNSTR RGTFSQELD AQLDDEDDN MHAVYPAVN RFSPSPRNSP RPSPKQSPRN SPRSRSPARG IEYSRVSPQP MISRLQQPRL SLQGHPTDLQ TSNVKNEEKL RRSLPNLSRT SNTQVDSVKS SRSDSNFQVP NGGIPRMQPQ ASAIPSPGKF RSPAAPSPLA LR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com