

Recombinant Human EHD3 Protein

Catalog #	Source	Reactivity	Applications
CST1499	E. coli	H	
Description	Recombinant Human EHD3 Protein is produced by E. coli expression system and the target gene encoding Leu291-Lys530 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	EHD3		
Alternative Names	EHD2; PAST3; EH domain-containing protein 3; PAST homolog 3		
Entrez Gene	30845		
SwissProt	Q9NZN3		
Purity			
Chemical Structure	LRKLNDLIKRLAKVHAYI ISSLKKEmps VFGKDNKKKE LVNNLAEIYG RIEREHQISP GDFPNLKRMQ DQLQAQDFSK FQPLKSKLLE VVDDMLAHDIAQLMVLVRQE ESQRPIQMVK GGAFEGTLHG PFGHGYGEA GEGIDDAEWV VARDKPMYDE IFYTLSPVDG KITGANAKKE MVRSKLPNSV LGKIWKLADI DKDGMLDDDE FALANHLLIKV KLEGHELPNE LP		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com