

Product Data Sheet

Recombinant Human CZIB Protein

Catalog #	Source	Reactivity	Applications
CST1501	E. coli	H	
Description	Recombinant Human CZIB Protein is produced by E. coli expression system and the target gene encoding Met1-Cys160 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	C1orf123		
Alternative Names	UPF0587 protein C1orf123		
Entrez Gene	54987		
SwissProt	Q9NWW4		
Purity			
Chemical Structure	MGKIALQLKA TLENITNLRP VGEDFRWYLK MKCGNCGEIS DKWQYIRLMD SVALKGGRGS ASMVQKCKLC ARENSIEILS STIKPYNAED NENFKTIVEF ECRGLEPVDF QPQAGFAAEG VESGTAFSDI NLQEKDWTDY DEKAQESVGI YEVTHQFVKC		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com