

Product Data Sheet

Recombinant Human TBC1D13 Protein

Catalog #	Source	Reactivity	Applications
CST1502	HEK293	H	
Description	Recombinant Human TBC1D13 Protein is produced by HEK293 expression system and the target gene encoding Pro81-Leu200 is expressed with C-terminus hFc & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TBC1D13		
Alternative Names	TBC1 domain family member 13		
Entrez Gene	54662		
SwissProt	Q9NVG8		
Purity			
Chemical Structure	PGIAKANMGV SREDVTFEDH PLNPNDSPRW NTYFKDNEVL LQIDKDVRRRL CPDISFFQRA TDYPCLLILD PQNEFETLRK RVEQTTLKSQ TVARNRSGVT NMSSPHKNSV PSSLNEYEVL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com