

Product Data Sheet

Recombinant Human NAA50 Protein

Catalog #	Source	Reactivity	Applications
CST1503	E. coli	H	
Description	Recombinant Human NAA50 Protein is produced by E. coli expression system and the target gene encoding Met1-Asn169 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	NAA50		
Alternative Names	MAK3; NAT13; NAT5; N-alpha-acetyltransferase 50; N-acetyltransferase 13; N-acetyltransferase 5; hNAT5; N-acetyltransferase san homolog; hSAN; NatE catalytic subunit		
Entrez Gene	80218		
SwissProt	Q9GZZ1		
Purity			
Chemical Structure	MKGSRIELGD VTPHNIKQLK RLNQVIFPVS YNDKFYKDVL EVGELAKLAY FNDIAVGAVC CRVDHSQNQK RLYIMTLGCL APYRRLGIGT KMLNHVLNIC EKDGTFDNIY LHVQISNESA IDFYRKFGFE IETKKNYYYK RIEPADAHVL QKNLKVPSGQ NADVQKTDN		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com