

Product Data Sheet

Recombinant Human CARD14 Protein

Catalog #	Source	Reactivity	Applications
CST1504	E. coli	H	
Description	Recombinant Human CARD14 Protein is produced by E. coli expression system and the target gene encoding Asp261-Ser480 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CARD14		
Alternative Names	CARMA2; Caspase recruitment domain-containing protein 14; CARD-containing MAGUK protein 2; Carma 2		
Entrez Gene	79092		
SwissProt	Q9BXL6		
Purity			
Chemical Structure	DEELNRLKEE NEKLRSITFS LAEKDILEQS LDEARGSRQE LVERIHSLRE RAVAAERQRE QYWEEKEQTL LQFQKSKMAC QLYREKVNAL QAQVCELQKE RDQAYSARDS ACREISQSLV EKDSLRRQVF ELTDQVCELR TQLRQLQAEP PGVLKQEART REPCPREKQR LVRMHAICPR DDSDCSLVSS TESQLLSDLS ATSSRELVDV FRSSSPAPPS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com