

Product Data Sheet

Recombinant Human DCUN1D1 Protein

Catalog #	Source	Reactivity	Applications
CST1507	E. coli	H	
Description	Recombinant Human DCUN1D1 Protein is produced by E. coli expression system and the target gene encoding Lys11-Gln250 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	DCUN1D1		
Alternative Names	DCUN1L1; RP42; SCCRO; DCN1-like protein 1; DCUN1 domain-containing protein 1; Defective in cullin neddylation protein 1-like protein 1; Squamous cell carcinoma-related oncogene		
Entrez Gene	54165		
SwissProt	Q96GG9		
Purity			
Chemical Structure	KVRQFMIFTQ SSEKTAVSCL SQNDWKLDVA TDNFFQNP EL YIRESVKGSL DRKKLEQLYN RYKDPQDENK IGIDGIQQFC DDLALDPASI SVLIIAWKFR AATQCEFSKQ EFMDGMTELG CDSIEKLKAQ IPKMEQELKE PGRFKDFYQF TFNFAKNPGQ KGLDLEMAIA YWNLVLNGRF KFLDLWNKFL LEHHKRSIPK DTWNLLLD FS TMIADDMSNY DEEGAWPVLI DD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com