

Product Data Sheet

Recombinant Human SH3KBP1 Protein

Catalog #	Source	Reactivity	Applications
CST1508	E. coli	H	
Description	Recombinant Human SH3KBP1 Protein is produced by E. coli expression system and the target gene encoding Glu331-Lys660 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SH3KBP1		
Alternative Names	CIN85; SH3 domain-containing kinase-binding protein 1; CD2-binding protein 3; CD2BP3; Cbl-interacting protein of 85 kDa; Human Src family kinase-binding protein 1; HSB-1		
Entrez Gene	30011		
SwissProt	Q96B97		
Purity			
Chemical Structure	EKEGNRPKKP PPPSAPVIKQ GAGTTERKHE IKKIPPERPE MLPNRTEEKE RPEREPKLDL QKPSVPAIPP KKPRPPKTNS LSRPGALPPR RPERPVGPLT HTRGDSPKID LAGSSLSGIL DKDLSDRSND IDLEGFDSVV SSTEKLSHPT TSRPKATGRR PPSQSLTSSS LSSPDIFDSP SPEEDKEEHI SLAHRGVDAS KKTSKTVTIS QVSDNKASLP PKPGTMAAGG GG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com