

Product Data Sheet

Recombinant Human LZIC Protein

Catalog #	Source	Reactivity	Applications
CST1511	E. coli	H	
Description	Recombinant Human LZIC Protein is produced by E. coli expression system and the target gene encoding Met1-Lys190 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LZIC		
Alternative Names	Protein LZIC; Leucine zipper and CTNNBIP1 domain-containing protein; Leucine zipper and ICAT homologous domain-containing protein		
Entrez Gene	84328		
SwissProt	Q8WZA0		
Purity			
Chemical Structure	MASRGKTETS KLKQNLEEQL DRLMQQLQDL EECREELDTD EYEETKKETL EQLSEFNDSL KKIMSGNMTL VDELSGMQLA IQAAISQAFK TPEVIRLFAK KQPGQLRTRL AEMDRDLMVG KLERDLYTQQ KVEILTALRK LGEKLTADDE AFLSANAGAI LSQFEKVSTD LGSGDKILAL ASFEVEKTKK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com