

Product Data Sheet

Recombinant Human GIPC3 Protein

Catalog #	Source	Reactivity	Applications
CST1512	E. coli	H	
Description	Recombinant Human GIPC3 Protein is produced by E. coli expression system and the target gene encoding Met1-Gly312 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	GIPC3		
Alternative Names	C19orf64; PDZ domain-containing protein GIPC3		
Entrez Gene	126326		
SwissProt	Q8TF64		
Purity			
Chemical Structure	MEGAAAREAR GTETPRASAP PPAPSEPPAA PRARPRLVFR TQLAHGSPTG KIEGFTNVRE LYAKIAEAFG IAPTEILFCT LNSHKVDMQK LLGGQIGLED FIFAHVRGET KEVEVTKTED ALGLTITDNG AGYAFIKRIK EGSIIINRIEA VCVGDSIEAI NDHSIVGCRH YEVAKMLREL PKSQPFTLRL VQPKRAFDMI GQRSRSSKCP VEAKVTSGRE TLRLRSGGAA TV		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com