

Product Data Sheet

Recombinant Human MRAP Protein

Catalog #	Source	Reactivity	Applications
CST1514	E. coli	H	
Description	Recombinant Human MRAP Protein is produced by E. coli expression system and the target gene encoding Ser81-Gly150 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	MRAP		
Alternative Names	C21orf61; FALP; Melanocortin-2 receptor accessory protein; B27; Fat cell-specific low molecular weight protein; Fat tissue-specific low MW protein		
Entrez Gene	56246		
SwissProt	Q8TCY5		
Purity			
Chemical Structure	SHGLNLHLCI QKCLPCHREP LATSQAQASS VEPGSRTGPD QPLRQESSST LPLGGFQTHP TLLWELTLNG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com