

Product Data Sheet

Recombinant Human TOR1AIP2 Protein

Catalog #	Source	Reactivity	Applications
CST1515	E. coli	H	
Description	Recombinant Human TOR1AIP2 Protein is produced by E. coli expression system and the target gene encoding Glu11-Ile210 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TOR1AIP2		
Alternative Names	IFRG15; LULL1; Torsin-1A-interacting protein 2; Luminal domain-like LAP1		
Entrez Gene			
SwissProt	Q8NFAQ		
Purity			
Chemical Structure	EDSQKDLEND PSVNSQAQET TIIASNAEEA EILHSACGLS KDHQEVETEG PESADTGDKS ESPDEANVGK HPKDKTEDEN KQSFLDGGKG HHLPSENLGK EPLDPDPSHS PSDKVGRADA HLGSSSVALP KEASDGTGAS QEPPTTDSQE AQSPGHSSAG QEGEDTLRRR LLAPEAGSHP QQTQKLEEI ENAQDTMRQI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com