

Recombinant Human ARFGAP2 Protein

Catalog #	Source	Reactivity	Applications
CST1516	E. coli	H	
Description	Recombinant Human ARFGAP2 Protein is produced by E. coli expression system and the target gene encoding Leu231-Glu390 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ARFGAP2		
Alternative Names	ZNF289; ADP-ribosylation factor GTPase-activating protein 2; ARF GAP 2; GTPase-activating protein ZNF289; Zinc finger protein 289		
Entrez Gene	84364		
SwissProt	Q8N6H7		
Purity			
Chemical Structure	LGAQKVSSQS FSEIERQAQV AEKLREQQAA DAKKQAEESM VASMRLAYQE LQIDRKKEEK KLQNLEGKKR EQAERLGMGL VSRSSVSHSV LSEMQUIEQE TPVSAKSSRS QLDFDDVGT FASGPPKYKD NPFSLGESFG SRWDTDAAWG MDRVEEKEPE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com