

Product Data Sheet

Recombinant Human CRLF3 Protein

Catalog #	Source	Reactivity	Applications
CST1519	E. coli	H	
Description	Recombinant Human CRLF3 Protein is produced by E. coli expression system and the target gene encoding Met1-Phe442 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	CRLF3		
Alternative Names	CREME9; CRLM9; CYTOR4; P48; Cytokine receptor-like factor 3; Cytokine receptor-like molecule 9; CREME-9; Cytokine receptor-related protein 4; Type I cytokine receptor-like factor p48		
Entrez Gene	51379		
SwissProt	Q8IUI8		
Purity			
Chemical Structure	MRGAMELEPE LLLQEARENV EAAQSYRREL GHRLEGLREA RRQIKESASQ TRDVLKQHFN DLKGTLGKLL DERLVTLTLLQE VDTIEQETIK PLDDCQKLIE HGVNTAEDLV REGEIAMLGG VGEENEKLWS FTKKASHIQL DSLPEVPLLVDVPCLSAQLD DSILNIVKDH IFKHGTVASR PPVQIEELIE KPGGIIVRWC KVDDDFTAQD YRLQFRKCTS NHFEDVYVGS ET		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com