

Recombinant Human HUWE1 Protein

Catalog #	Source	Reactivity	Applications
CST1520	E. coli	H	
Description	Recombinant Human HUWE1 Protein is produced by E. coli expression system and the target gene encoding Phe4001-Gly4370 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HUWE1		
Alternative Names	KIAA0312; KIAA1578; UREB1; E3 ubiquitin-protein ligase HUWE1; ARF-binding protein 1; ARF-BP1; HECT, UBA and WWE domain-containing protein 1; Homologous to E6AP carboxyl terminus homologous protein 9; HectH9; Large structure of UREB1; LASU1; Mcl-1 ubiquiti		
Entrez Gene	10075		
SwissProt	Q7Z6Z7		
Purity			
Chemical Structure	FRQELERLDE GLRKEDMAVH VRRDHVFEDS YRELHRKSPE EMKNRLYIVF EGEEGQDAGG LLREWYMIIS REMFNPMYAL FRTSPGDRVT YTINPSSHCHN PNHLSYFKFV GRIVAKAVYD NRLLECIFYTR SFYKHILGKS VRYTDMESD YHFYQGLVYL LENDVSTLGY DLTFTSTEVQE FGVCEVRDLK PNGANILVTE ENKKEYVHLV CQMRMTGAIR KQLAAFLEGF YE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com