

Product Data Sheet

Recombinant Human CRACD Protein

Catalog #	Source	Reactivity	Applications
CST1522	E. coli	H	
Description	Recombinant Human CRACD Protein is produced by E. coli expression system and the target gene encoding Pro801-Thr1180 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIAA1211		
Alternative Names	Uncharacterized protein KIAA1211		
Entrez Gene	57482		
SwissProt	Q6ZU35		
Purity			
Chemical Structure	PSLPYPPQKV VAHTEFTTSS DSETANGIAK PDPVMPGGEE KASPFGIKLR RTNYSLRFNC DQQAEEQKKKK RHSSTGDSAD AGPPAAGSAR GEKEMEGVAL KHGPSLPQER KQAPSTRRDS AEPSSSRVSP VAHPGPPPAS SQTPAPEHDK AANKMPLAQK PALAPKPTSQ TPPASPLSKL SRPYLVLLS RRAGRPDPEP SEPSKEDQES SDRRPPSPPG PEERKGQKRD EE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com