

Product Data Sheet

Recombinant Human SKIC3 Protein

Catalog #	Source	Reactivity	Applications
CST1524	E. coli	H	
Description	Recombinant Human SKIC3 Protein is produced by E. coli expression system and the target gene encoding Lys11-Ile230 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TTC37		
Alternative Names	KIAA0372; Tetratricopeptide repeat protein 37; TPR repeat protein 37; SKI3 homolog; Ski3; Tricho-hepatic-enteric syndrome protein; Thespin		
Entrez Gene	9652		
SwissProt	Q6PGP7		
Purity			
Chemical Structure	KSARDAIRNK EYKEALKHCK TVLKQEKNNY NAWVFIGVAA AELEQPDQAAQ SAYKKAAELE PDQLLAWQGL ANLYEKYNHI NAKDDLPGVY QKLLDLYESV DKQKWCDVCK KLVDLYYQEK KHLEVARTWH KLIKTRQEQG AENEELHQLW RKLTQFLAES TEDQNNETQQ LLFTAFENAL GLSDKIPSED HQVLYRHFIQ SLSKFPHESA RLKKACEGMI		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com