

Product Data Sheet

Recombinant Human LCE3C Protein

Catalog #	Source	Reactivity	Applications
CST1525	HEK293	H	
Description	Recombinant Human LCE3C Protein is produced by HEK293 expression system and the target gene encoding Gln11-Ser90 is expressed with N-terminus GST & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LCE3C		
Alternative Names	LEP15; SPRL3A; Late cornified envelope protein 3C; Late envelope protein 15; Small proline-rich-like epidermal differentiation complex protein 3A		
Entrez Gene	353144		
SwissProt	Q5T5A8		
Purity			
Chemical Structure	QPPPSCPSPK CPPKSPAQCL PPPSSDCALS SGGCGPSSES GCCLSHHRHF RSHQCRRQRS NSCDRGSGQQ GGGSCRGHGS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com