

Product Data Sheet

Recombinant Human ADIRF Protein

Catalog #	Source	Reactivity	Applications
CST1526	E. coli	H	
Description	Recombinant Human ADIRF Protein is produced by E. coli expression system and the target gene encoding Met1-Lys76 is expressed with N-terminus SUMO & C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ADIRF		
Alternative Names	AFRO; APM2; C10orf116; Adipogenesis regulatory factor; Adipogenesis factor rich in obesity; Adipose most abundant gene transcript 2 protein; Adipose-specific protein 2; apM-2		
Entrez Gene	10974		
SwissProt	Q15847		
Purity			
Chemical Structure	MASKGLQDLK QQVEGTAQEA VSAAGAAAQQ VVDQATEAGQ KAMDQLAKTT QETIDKTANQ ASDTFSGIGK KFGLLK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com