

Product Data Sheet

Recombinant Human SELENBP1 Protein

Catalog #	Source	Reactivity	Applications
CST1531	E. coli	H	
Description	Recombinant Human SELENBP1 Protein is produced by E. coli expression system and the target gene encoding Met1-Ile472 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	SELENBP1		
Alternative Names	SBP; Selenium-binding protein 1; 56 kDa selenium-binding protein; SBP56; SP56		
Entrez Gene	8991		
SwissProt	Q13228		
Purity			
Chemical Structure	MATKCGNCGP GYSTPLEAMK GPREEIVYLP CIYRNTGTEA PDYLATVDVD PKSPQYCQVI HRLPMPNLKD ELHHSWNTC SSCFGDSTKS RTKLVLSLI SSRIYVVDVG SEPRAPKLHK VIEPKDIHAK CELAFLHTSH CLASGEVMIS SLGDVKGNKG GGFVLLDGET FEVKGTWERP GGAAPLG YDF WYQPRHN VMI STEWAAP NVL RDGFNPADVE AGLYGSHLYV WD		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com