

Product Data Sheet

Recombinant Human FLII Protein

Catalog #	Source	Reactivity	Applications
CST1532	E. coli	H	
Description	Recombinant Human FLII Protein is produced by E. coli expression system and the target gene encoding Trp941-Ala1260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	FLII		
Alternative Names	FLIL; Protein flightless-1 homolog		
Entrez Gene	2314		
SwissProt	Q13045		
Purity			
Chemical Structure	WVPVEYEEEE KKEDKEEKAE GKEGEEATAE AEEKQPEEDF QCIVYFWQGR EASNMGWLTF TFSLQKKFES LFP GKLEVVR MTQQQENPKF LSHFKRKFI HRGKRKAVQG AQQPSLYQIR TNGSALCTRC IQINTDSSL NSEFCILKV PFESEDNQGI VYAWVGRASD PDEAKLAEDI LNTMFDTSYS KQVINEGEEP ENFFWVGIGA QKPYDDDAEY MKHTRLFRCS NE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com