

Product Data Sheet

Recombinant Human hnRNP U Protein

Catalog #	Source	Reactivity	Applications
CST1535	E. coli	H	
Description	Recombinant Human hnRNP U Protein is produced by E. coli expression system and the target gene encoding Leu11-Tyr260 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	HNRNPU		
Alternative Names	HNRPU; SAFA; U21.1; Heterogeneous nuclear ribonucleoprotein U; hnRNP U; Scaffold attachment factor A; SAF-A; p120; pp120		
Entrez Gene	3192		
SwissProt	Q00839		
Purity			
Chemical Structure	LKVSELKEEL KKRRLSDKGL KAELMERLQA ALDDEEAGGR PAMEPGNGSL DLGGDSAGRS GAGLEQEAAA GGDEEEEEEE EEEEGISALD GDQMELGEEN GAAGAADSGP MEEEEAASED ENGDDQGFQE GEDELGDEEE GAGDENGHGE QQPQPPATQQ QQPQQRGAA KEAAGKSSGP TSLFAVTVP PGARQGQQQA GGKKKAEGGG GGGRPGAPAA GDGKTEQKGG DK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com