

Recombinant Human UBASH3A Protein

Catalog #	Source	Reactivity	Applications
CST1537	E. coli	H	
Description	Recombinant Human UBASH3A Protein is produced by E. coli expression system and the target gene encoding Val401-Ala650 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	UBASH3A		
Alternative Names	STS2; Ubiquitin-associated and SH3 domain-containing protein A; Cbl-interacting protein 4; CLIP4; Suppressor of T-cell receptor signaling 2; STS-2; T-cell ubiquitin ligand 1; TULA-1		
Entrez Gene	53347		
SwissProt	P57075		
Purity			
Chemical Structure	VRHGERVDQI FGKAWLQQCS TPDGKYRPPD LNFPCSLPRR SRGIKDFEND PPLSSCGIFQ SRIAGDALLD SGIRISSVFA SPALRCVQTA KLILEELKLE KIKIRVEPG IFEWTKWEAG KTTPTLMSLE ELKEANFNID TDYRPAFPLS ALMPAESYQE YMDRCTASMV QIVNTCPQDT GVILIVSHGS TLDSCSTRPLL GLPPRECGDF AQLVRKIPSL GCMCFCEENKE EG		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com