

Recombinant Human LRPPRC Protein

Catalog #	Source	Reactivity	Applications
CST1545	E. coli	H	
Description	Recombinant Human LRPPRC Protein is produced by E. coli expression system and the target gene encoding Ile61-Tyr320 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	LRPPRC		
Alternative Names	LRP130; Leucine-rich PPR motif-containing protein, mitochondrial; 130 kDa leucine-rich protein; LRP 130; GP130		
Entrez Gene	10128		
SwissProt	P42704		
Purity			
Chemical Structure	IAAKEKDIQE ESTFSSRKIS NQFDWALMRL DLSVRRTGRI PKLLQKVFN DTCRSGGLGG SHALLLRSC GSLLPELKLE ERTEFAHRIW DTLQKLGAVY DVSHYNALLK VYLQNEYKFS PTDFLAKMEE ANIQPNRVTY QRLIASYCNV GDIEGASKIL GFMKTKDLPV TEAVFSALVT GHARAGDMEN AENILTVMRD AGIEPGPDY LALLNAYA EK GDIDHVKQTL EK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com