

Recombinant Human TEC Protein

Catalog #	Source	Reactivity	Applications
CST1546	E. coli	H	
Description	Recombinant Human TEC Protein is produced by E. coli expression system and the target gene encoding Leu11-Arg340 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TEC		
Alternative Names	PSCTK4; Tyrosine-protein kinase Tec		
Entrez Gene	7006		
SwissProt	P42680		
Purity			
Chemical Structure	LIKRSQQKKK TSPLNYKERL FVLTKSMLTY YEGRAEKKYR KGFIDVSKIK CVEIVKNDDG VIPCQNKYPF QVVHDANTLY IFAPSPQSRD LWVKKLKEEI KNNNNIMIKY HPKFWTDGSY QCCRQTEKLA PGCEKYNLFE SSIRKALPPA PETKKRRPPP PIPEEEDNS EEIVVAMYDF QAAEGHDLRL ERGQEYLILE KNDVHWWRAR DKYGNEGYIP SNYVTGKKS NL		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com