

Recombinant Human ELK3 Protein

Catalog #	Source	Reactivity	Applications
CST1547	E. coli	H	
Description	Recombinant Human ELK3 Protein is produced by E. coli expression system and the target gene encoding Leu11-Lys290 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	ELK3		
Alternative Names	NET; SAP2; ETS domain-containing protein Elk-3; ETS-related protein ERP; ETS-related protein NET; Serum response factor accessory protein 2; SAP-2; SRF accessory protein 2		
Entrez Gene	2004		
SwissProt	P41970		
Purity			
Chemical Structure	LLQLLLDQKH EHLICWTSND GEFKLLKAEV VAKLWGLRKN KTNMNYDKLS RALRYYYDKN IIKKVIGQKF VYKFVSFPEI LKMDPHAVEI SRESLLQDS DCKASPEGRE AHKHGLAALR STSRNEYIHS GLYSSFTINS LQNPPDAFKA IKTEKLEPP EDSPPVVEVR TVIRFVTNKT DKHVTRPVVS LPSTSEAAAA SAFLASSVSA KISSLMPLNA ASISSASPFS SR		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com