

Product Data Sheet

Recombinant Human TACC2 Protein

Catalog #	Source	Reactivity	Applications
CST1565	E. coli	H	
Description	Recombinant Human TACC2 Protein is produced by E. coli expression system and the target gene encoding Leu2751-Glu2940 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TACC2		
Alternative Names	Transforming acidic coiled-coil-containing protein 2; Anti-Zuai-1; AZU-1		
Entrez Gene	10579		
SwissProt	O95359		
Purity			
Chemical Structure	LQIARAEIIT KEREVSEWKD KYEESRREVM EMRKIVAEYE KTIAQMIEDE QREKSVSHQT VQQLVLEKEQ ALADLNSVEK SLADLFRRYE KMKEVLEGFR KNEEVLRCA QEYLSRVKKE EQRYQALKVH AEEKLDRANA EIAQVRGKAQ QEQAHHQASL RKEQLRVDAL ERTLEQKNKE IEELTKICDE		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com