

Product Data Sheet

Recombinant Human KIAA0513 Protein

Catalog #	Source	Reactivity	Applications
CST1569	E. coli	H	
Description	Recombinant Human KIAA0513 Protein is produced by E. coli expression system and the target gene encoding Met1-Glu411 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	KIAA0513		
Alternative Names	Uncharacterized protein KIAA0513		
Entrez Gene	9764		
SwissProt	O60268		
Purity			
Chemical Structure	METPEVPVGS LIDFGPEAPT SSPLEAPPPV LQDGDGSLGD GASESETTES ADSENDMGES PSHPSWDQDR RSSSNESFSS NQSTESTQDE ETLALRDFMR GYVEKIFSGG EDLDQEEKAK FGEYCSENG KGREWFARYV SAQRCNSKCV SEATFYRLVQ SFAVVLFECH QMDDFGPAKN LMTMCFTYYH IGKPQLLPPE SREKPAGSID SYLKSANSWL AEKKDIAERL LK		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com