

Product Data Sheet

Recombinant Human TRAF3IP2 Protein

Catalog #	Source	Reactivity	Applications
CST1570	E. coli	H	
Description	Recombinant Human TRAF3IP2 Protein is produced by E. coli expression system and the target gene encoding Asn11-Leu400 is expressed with C-terminus His tag.		
Form	Liquid in a 0.2 µm filtered solution of PBS, pH 7.4.		
Gene Symbol	TRAF3IP2		
Alternative Names	C6orf2; C6orf4; C6orf5; C6orf6; Adapter protein CIKS; Connection to IKK and SAPK/JNK; Nuclear factor NF-kappa-B activator 1; ACT1; TRAF3-interacting protein 2		
Entrez Gene	10758		
SwissProt	O43734		
Purity			
Chemical Structure	NRSIPVEVDE SEPYPSQLLK PIPEYSPEEE SEPPAPNIRN MAPNSLSAPT MLHNSSGDFS QAHSTLKLAN HQRPVSRQVT CLRTQVLEDS EDSFCRRHPG LGKAFPSGCS AVSEPASESV VGALPAEHQF SFMEKRNQWL VSQLSAASPD TGHDSKSDQ SLPNASADSL GGSQEMVQRP QPHRNAGLD LPTIDTGYDS QPQDVLGIRQ LERPLPLTSV CYPQDLPRPL RS		
Quality Control			
Directions for Use	Always centrifuge tubes before opening. Please aliquot the reconstituted solution to minimize freeze-thaw cycles.		
Storage/Stability	Store it at -20°C to -80°C for one year. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.		

Application key: E- ELISA, WB- Western blot, IH- Immunohistochemistry, IF- Immunofluorescence, FC- Flow cytometry, IC- Immunocytochemistry, IP- Immunoprecipitation, ChIP- Chromatin Immunoprecipitation, EMSA- Electrophoretic Mobility Shift Assay, BL- Blocking, SE- Sandwich ELISA, CBE- Cell-based ELISA, RNAi- RNA interference

Species reactivity key: H- Human, M- Mouse, R- Rat, B- Bovine, C- Chicken, D- Dog, G- Goat, Mk- Monkey, P- Pig, Rb- Rabbit, S- Sheep, Z- Zebrafish

COHESION BIOSCIENCES LIMITED

WEB
www.cohesionbio.com

ORDER
order@cohesionbio.com

SUPPORT
techsupport@cohesionbio.com

CUSTOM
custom@cohesionbio.com